

Infection Prevention and Control

Mechanical Lift Disinfecting Audit Tool

Date/Time of Observation						
Initials of Observer						
Name/Title of Caregiver						
Staff state reason(s) for disinfection of mechanical lifts. Yes or No						
Lift is disinfected prior to use according to the manufacturer's guidelines, policies and procedures (P/P). Yes or No						
Staff perform hand hygiene (HH) and wear gloves prior to disinfecting (follow directions for use of disinfectant). Yes or No						

Appropriate disinfectant is used. Yes or No						
Staff clean boom, sling bar, control, handles (touch surfaces). Yes or No						
Staff demonstrate waiting appropriate dry time before use (2 minutes). Yes or No						
Staff disinfect lift after each resident use and when obtaining a lift from storage (HH, gloves). Yes or No						
Staff state process for refill of disinfectant wipes (where extra wipes are located, replacing if empty). Yes or No						
Notes:						